



NO.1 CHOICE OF HOSPITALS



Breastfeeding Products

MODELLED ON NATURE

Overview

Breastfeeding – naturally good	4
Breastfeeding Products	5
Breastpump Overview	6
2-Phase Expression®	7
Breastpumps of the latest generation	8 / 9
Classic Breastpumps	10 / 11
Storing and transporting breastmilk	12
Breastshields	13
Nipple Care	14
Milk Leakage	15
Special Feeding Devices	16 / 17
Medela – «Made in Switzerland»	18

Breastfeeding – naturally good

MOTHER'S MILK IS BEST

Breastfeeding is one of the most beautiful and natural experiences a mother can share with her baby. A unique time when they feel intensely close and intimate. And both are completely contented. Breastfeeding also benefits the health of both mother and baby in many different ways.



Ideal for mother and baby

With her breastmilk the mother gives her baby all the important nutrients and antibodies it needs for its development. Breastfeeding protects both mother and baby from illnesses and promotes lifelong good health. Breastfeeding is also practical. The milk is always at the right temperature and immediately available whenever the baby gets hungry.

Recommended by WHO and Unicef

Nature knows what is good. Which is why the World Health Organization (WHO), Unicef, lactation consultants and paediatricians recommend that for the first six months babies should only be fed on breastmilk. And that partial breastfeeding should ideally continue into their second year.

Breastfeeding is good for mother and baby – easy and natural.

Breastfeeding Products

MODELLED ON NATURE

Mothers sometimes need a little gentle support with breastfeeding. Medela knows their individual needs and gives them its total attention. That is why mothers, lactation consultants, other specialists as well as distinguished research teams are involved when its innovative breastfeeding products are developed.

Researching nature

Medela is strongly dedicated to research. For further development of its breastfeeding products it collaborates with research teams at the University of Western Australia and other selected teams elsewhere. That way, latest scientific knowledge is included in the product development process to benefit the mother and her baby.

Gentle support

Breastpumps can be used to support breastfeeding in special situations. For example, when a premature baby is still in hospital, or the mother does not produce enough milk, or too much, or has sore nipples.

But they are also part of a modern lifestyle. You want to breastfeed for as long as possible but still be independent or go back to work? Mothers can express their milk at home or when they are out and about. And the baby can enjoy its mother's milk while its father cuddles it at home.

That way, mothers and babies can enjoy breastfeeding for a long time.



Breastpump Overview

INDIVIDUAL SOLUTIONS TO SUIT MOTHER AND BABY

One mother wants a discreet breastpump she can take with her when she goes out. Another wants to be able to express both breasts simultaneously. Medela has the wide range of products it takes to meet these individual needs. Medela offers innovative breastpumps for every lifestyle and every budget.

Choosing the right breastpump

	For hospitals and rental electric		For personal use electric				For personal use manual	
	Double Pumping		Single Pumping					
	Symphony	Lactina	Pump In Style Advanced	Mini Electric Plus	Swing	Mini Electric	Harmony	Manual Breastpump
Occasional Use								
A missed feeding	●	●	●	●	●	●	●	●
An evening out from baby	●	●	●	●	●	●	●	●
Working part-time	●	●	●	●	●	●	●	●
Daily Use								
Working full-time	●	●	●		●			
Medical Reasons								
Premature or hospitalized baby	●	●						
Low milk supply	●	●						
Latch-on problems or breast infection	●	●	●		●			
Sore nipples or engorgement	●	●	●	●	●	●		
Drawing out flat or inverted nipples	●	●	●	●	●	●	●	●

1st choice ● 2nd choice ●

Medela Breastpumps are recommended by lactation consultants and midwives all over the world and support extended breastfeeding.



2-Phase Expression®

THE GOLD STANDARD – EXCLUSIVELY FROM MEDELA


Breastmilk is best for the baby. And the baby's natural sucking rhythm is best for the mother. In intensive pioneer work, Medela has converted the baby's intuitive knowledge into technological know-how. The result is 2-Phase Expression® – the new standard for unsurpassed naturalness, comfort and efficiency in breastpumps.

Closer to nature than ever

Babies change the way they suckle at a mother's breast during breastfeeding. Initially, they suck rapidly to stimulate the milk ejection reflex (MER) and induce milk «let-down». Once the milk starts flowing, they settle into a more relaxed sucking speed. The Medela 2-Phase Expression simulates this natural nursing rhythm more precisely than any other system. It makes pumping especially gentle, comfortable and efficient.

Only the baby is better!

2-Phase Expression® – The Original



«2-Phase Expression is a true innovation. Solidly researched, perfectly implemented.»

Ines F., Lactation Consultant



Breastpumps of the latest generation

NATURE AND TECHNOLOGY PERFECTLY COMBINED

Medela makes 2-Phase Expression® the standard. Initially only available with the Symphony, this innovation is now presented in three further models. For every mother and every situation.

Symphony®

This breastpump for hospitals and home rental gives the mother the most natural of feelings along with unique technical refinement. Ideal for long-term and frequent pumping needs.

- Single or double pumping
- Whisper quiet
- Overflow protection
- Mobile stand available



Pump In Style® Advanced

With this electric breastpump in a trendy backpack the active and working mother is fully equipped for pumping – any time, any place, every day. Ideal for frequent and extended use.

- Single or double pumping
- Removable motor unit and Cooler Carrier
- Electric or battery-powered



Swing™

An electric breastpump which perfectly adapts to the daily routine of a dynamic mother. Ideal for occasional or frequent pumping – handy, trendy and quiet.

- Ideal for out and about
- Gently massaging SoftFit™ breastshield
- Electric or battery-powered
- Cooler Carrier available



Harmony™

This manual breastpump is ideal for mothers who only express milk occasionally and do not want to miss the usual feel of their baby.

- Ergonomic swivel handle
- Gently massaging SoftFit™ breastshield
- Fewer parts than other manual breastpumps
- Cooler Carrier available



Classic Breastpumps

PROVEN BY MOTHERS WORLDWIDE

Mothers have been relying on Medela breastfeeding products for many years. Medela's classic breastpumps also meet their countless individual needs. For every situation and every budget.

Lactina™ Electric Plus

Lactina™ Select

These electric breastpumps are ideal for long-term and regular use in hospitals or as rented pumps for home use.

- Single or double pumping
- The accessory set can also function as a manual breastpump
- Mobil stand available

Lactina Select

On the Lactina Select both sucking strength and pumping speed are adjustable



Lactina Electric Plus



Mini Electric™ Plus

This fully automatic breastpump allows double pumping with all its advantages. Ideal for expressing breastmilk occasionally when at home or on the go.

- Single or double pumping
- Electric or battery-powered
- Cooler Carrier available

Mini Electric™

This small electric breastpump offers optimal comfort for expressing breastmilk at home or on the go. Ideal for short-term and occasional use.

- Easy to clean and assemble
- Electric or battery-powered
- Cooler Carrier available

Millions already sold



Manual Breastpump

The Manual Breastpump offers convenient manual pumping. Ideal for occasional use.

- Easy to clean and assemble
- Mini Electric Conversion Kit also available
- Cooler Carrier available

Storing and transporting breastmilk

CONVENIENT AND EASY

Medela offers breastfeeding mothers a variety of accessories that make pumping, collecting and storing breastmilk easier and more comfortable.



Breastmilk Bottles

Collect, store and freeze breastmilk.

- Milk may be used directly from bottle when thawed
- Can be used with all Medela breastpumps
- Made of polypropylene, the safest plastic for storing breastmilk

Breastmilk Bags

Collect, store and freeze breastmilk.

- Pre-sterilized for maximum hygiene
- Convenient easy pumping directly into bag
- Built-in, easy-to-open pouring spout
- Safe and hygienic storage at home or in hospitals



Cooler Carrier

To cool, store and transport breastmilk.

- Offers additional room for the breastpumps Swing, Harmony, Mini Electric Plus, Mini Electric or Manual Breastpump
- Removable cooler bag

Breastshields

THE RIGHT FIT FOR EVERY MOTHER

Every mother is unique – one breastshield does not fit all. Medela offers a variety of breastshields to meet each mother's needs.

PersonalFit™ Breastshield

Customized two-piece system to accommodate the needs of more mothers.

- For larger than average nipples
- Available in two tunnel sizes (27 and 30 mm)



SoftFit™ Breastshield

The SoftFit Breastshield is soft, flexible and comfortably massages the breast and areola.

- Imitates a baby's sucking action
- For maximum comfort

Comfort Breastshield

For higher comfort while pumping.

- Imitates a baby's sucking action
- Made of soft silicone



Nipple Care

EASY TO USE – EFFECTIVE SUPPORT

Sore nipples are one of the main reasons for breastfeeding difficulties. These products can provide relief and support. Mothers with flat or inverted nipples can breastfeed, too – Nipple Formers help prepare their nipples for breastfeeding.

PureLan™ 100

Soothes sensitive or dry nipples, dry baby skin and dry skin.

- All-natural, hypoallergenic lanolin
- Replenishes the natural oils depleted by frequent sucking
- No need to remove before breastfeeding



Nipple Formers

A non-invasive method to prepare flat or inverted nipples for breastfeeding.

- Soft silicone backs for a snug, comfortable fit
- Vented fronts to allow air circulation
- Contoured for less visibility through clothing



Breastshells

Protect sore or cracked nipples.

- Soft silicone backs for a snug, comfortable fit
- Vented fronts to allow air circulation
- Contoured for less visibility through clothing



Contact™ Nipple Shields

For temporary use when breastfeeding with sore, cracked or flat nipples. Help with latch-on difficulties.

- Open section for baby to stay in contact with mother's familiar smell
- Made of thin, soft silicone
- Neutral in taste
- Available in two sizes



Milk Leakage

HYGIENIC AND DISCREET COLLECTION AND ABSORPTION

For breastmilk that leaks while breastfeeding or between feeds, Medela offers a range of helpful products.

Disposable Bra Pads

Absorb milk leakage.

- Unique moisture-absorbing materials
- Adhesive tape keeps pad in place
- Contoured for a natural, smooth fit



Washable Bra Pads

Absorb milk leakage.

- Allow air circulation around the nipple
- Contoured shape for a comfortable fit
- Washable and reusable



Milk Collection Shells

Collect leaking milk.

- Soft silicone backs for a snug, comfortable fit
- Easy emptying through convenient spout
- Contoured for less visibility through clothing



Special Feeding Devices

PROFESSIONAL SUPPORT FOR SPECIAL NEEDS

Although most mothers and babies can breastfeed without special assistance, sometimes occasional help is needed. Premature babies, failure-to-thrive babies and babies with specific impairments, such as neurological challenges, may require extra attention. Infants with severe feeding problems can also benefit from Medela's specialized feeding products.

The products described here are designed to help overcome circumstances that might otherwise impede a successful breastfeeding experience. Seeking consultation with a lactation professional is recommended when using these products.

SoftCup™ Advanced Cup Feeder



Eases transition to breastfeeding. Medela's SoftCup Feeder is designed for use when health care professionals recommend cup feeding or for feeding breastfed babies for whom nipple confusion is a concern.

- Control valve and self-filling reservoir for more feeding control
- Special design prevents spillage of breastmilk
- Contoured to fit even the smallest mouth
- Made of soft, pliable silicone
- No risk of nipple confusion

Disposable Baby Cup Feeder



One method of feeding babies with breastmilk when nursing is not possible.

- Feeding lip for greater control while feeding
- Graduated markings



Please seek consultation when using these products.



Haberman Feeder

Helps ease feeding complications. For babies with nursing difficulties, particularly children with cleft lip/palate, weak babies or babies with Down's syndrome or other neurological difficulties.

- Variable flow rate to suit the baby's needs
- Sensitive to the weakest feeding effort
- One-way valve prevents flooding or overwhelming the baby

FingerFeeder



To teach babies proper suckling. Can also be used for feeding small amounts of breastmilk. The soft silicone tip is attached to a syringe.

Supplemental Nursing System (SNS)



For long-term supplemental feedings at the breast. Keeps the baby on the breast, ensuring the stimulation to maintain lactation and the intimacy between mother and baby.

- Adjustable flow rate system to meet individual needs
- Colour-coded tubes for different feeding rates
- Graduated reservoir bottle



Please seek consultation when using these products.

Medela – «Made in Switzerland»

A HIGH QUALITY STANDARD

A high quality standard and the best possible safety for users are integral components of Medela's corporate philosophy. We develop and manufacture our products with a certified quality management system in accordance with ISO 9001 and ISO 13485 as well as with the American FDA guidelines for GMP (Good Manufacturing Practices). Our products bear the CE mark in accordance with the Council Directive 93/42/EEC for medical products.



«We strive to provide the best products for mothers and babies. For us, quality does not stop at fulfilling the highest quality standards but also means steadily improving the breastfeeding mother's quality of life. That means constantly increasing her mobility and flexibility so she can always give her baby the best that nature can offer – her own milk.»

Michael Larsson, Chairman

Canada

Medela, Inc.
4090B Sladeview Crescent, Unit 2
Mississauga, Ontario
Canada, L5L 5Y5
Phone +1 905 608 7272
Fax +1 905 608 8720
e-mail customer.service@medela.ca
www.medela.com

France

Medela France Sarl
29, avenue des Grenots
91150 Etampes
Phone +33 (0)1 69 16 10 30
Fax +33 (0)1 69 16 10 32
e-mail info@medela.fr
www.medela.fr

Germany

Medela Medizintechnik
GmbH & Co. Handels KG
Postfach 1148
85378 Eching
Phone +49 (0)89 31 97 59-0
Fax +49 (0)89 31 97 59 99
e-mail info@medela.de
www.medela.de

Italy

Medela Italia Srl
Via Stelloni 39/O
40010 Sala Bolognese
Phone +39 051 681 52 76
Fax +39 051 681 52 77
e-mail info@medela.it
www.medela.it

Japan

Medela K.K.
Fuji Chuo Building 5F
2-1-7, Shintomi
Chuo-ku, Tokyo 104-0041
Phone +81 3 3553 3571
Fax +81 3 3553 3573
e-mail info@medela.co.jp
www.medela.co.jp

Netherlands

Medela Benelux B.V.
Dintel 4
5422 VT Gemert
Phone +31 0 492 36 24 23
Fax +31 0 492 36 72 80
e-mail info@medela.nl
www.medela.nl

Russia

OOO Medela
Fersmana Street 5A
117312 Moskau
Phone +7 (095) 980 6194
Fax +7 (095) 718 1700
e-mail info@medela-russia.ru
www.medela-russia.ru

Spain

Medela España
c/ Manuel Fernández Márquez, 49
08918 Badalona
Phone +34 93 320 59 69
Fax +34 93 320 55 31
e-mail info@medela.es
www.medela.es

Sweden

Medela Medical AB
Box 7266
187 14 Täby
Phone +46 (0)8 630 09 40
Fax +46 (0)8 630 09 48
e-mail info@medela.se
www.medela.se

United Kingdom

Medela UK Ltd.
Huntsman Drive
Northbank Industrial Park
Irlam, Manchester M44 5EG
Phone 0870 950 5994
Fax +44 (0)161 776 0444
e-mail info@medela.co.uk
www.medela.co.uk

USA

Medela, Inc.
P.O. Box 660
McHenry, IL. 60051-0660
Phone +1 815 363 1166
Fax +1 815 363 1246
e-mail info@medela.com
www.medela.com

Switzerland & International Sales

Medela AG, Medical Technology
Lättichstrasse 4b
6341 Baar/Switzerland
Phone +41 (0)41 769 51 51
Fax +41 (0)41 769 51 00
e-mail info@medela.ch
www.medela.ch